Novilase[®] Breast ILT Kit

Model #: 100010-01 and 100010-02

Instructions for Use

	<p>Novian Health Inc. 222 W. Merchandise Mart Plaza Suite 570 Chicago, IL 60654 USA 312-266-7200</p>
---	--

E.C. Rep Ltd

5 Fitzwilliam Square East
Dublin 2, D02 R744, Ireland
Tel: +353 89 225 1951
Email: info@ecrep.ie

UK Responsible Person

Responsible Person Ltd.
Summit House, 4-5 Mitchell Street
Edinburgh, EH6 7BD Scotland
+44 131 618 1115

info@responsible-person.co.uk

MedEnvoy Switzerland

Gotthardstrasse 28
6302 Zug
Switzerland

Table of Contents

1	Intended use of Novilase Therapy System	3
2	Indications for Use.....	3
3	Device Description & Sterilization	3
4	Contraindications.....	4
5	Warnings	4
6	Precautions	4
7	Residual Risks and undesirable side-effects.....	5
8	Storage and handling conditions.....	6
9	Novilase Procedure	6
10	Post-Novilase Procedure Instructions	9
11	Expected Clinical Benefit.....	9
12	Symbols	10

Before using the Novilase LTS 2.1 system with its accessories, you must read and fully understand the Novilase LTS 2.1 system User Manual. This document only provides instructions on using the Novilase laser and thermal probes (i.e., Novilase Breast Laser Therapy Kit) with the Novilase LTS 2.1 laser ablation system. Users must complete the user training course before using the system.

Novilase® is a Trademark of Novian Health Inc. in Europe, the United States, and other countries. Novian Health Inc. and other parties may also have trademark rights in other terms used herein.

1 Intended use of Novilase Therapy System

The Novilase® Laser Therapy System is intended for focal destruction of breast tumors by using laser-induced heat.

2 Indications for Use

Europe, UK, and Switzerland

The Novilase Laser Therapy System can be used for benign breast tumours up to 20 mm in the longest diameter.

The Novilase Laser Therapy System can be used for malignant breast tumours up to 20 mm in the longest diameter.

United States

The Novilase® Laser Therapy System (LTS-2.1) is indicated for the treatment of fibroadenomas of the breast, with tumour sizes up to 20 mm, and for general surgery procedures including incision, excision, and ablation of soft tissues; and coagulative necrosis and interstitial laser coagulation of soft tissues.

3 Device Description & Sterilization

The Novilase Laser Therapy Pack is a procedure pack containing the disposable components of the Novilase Laser Therapy System for use with the Novilase LTS 2.1 only and other CE marked supplies for the procedure. Included is:

1 sterile Novilase Breast ILT Kit containing:

- a) 1 Laser probe with 1 removable trocar stylet
- b) 1 Thermal probe
- c) 1 Laser fiber
- d) 1 Probe guide and 1 sterile BD leur-lok syringe, 30 mL for saline.

It is shipped with the probes, laser fiber, and probe guide sterilized with gamma radiation. This procedure pack is intended for single patient use only; discard all items after use.

Sites will be asked to supply sterile pressure monitoring tubing for the saline (72" male/female) and a sterile hemostasis valve (e.g., Merit Hemostasis High Pressure Y with Male Luer Lock, up to 6 FR).

The temperature sensors are capable of measuring temperature between 22 °C and 100°C with an accuracy of +/-1°C.

Model # 100010-01 contains 150 mm laser and thermal probes, and Model # 100010-02 contains 110 mm laser and thermal probes.

4 **Contraindications**

The Novilase Laser Therapy System is not indicated for:

- a) Pregnant or lactating women
- b) Women with uncorrectable coagulopathy
- c) A breast tumor that cannot be clearly visualized by ultrasound or MRI (if malignant)
- d) DCIS with microinvasion
- e) A cluster of microcalcifications that exceed 10 mm in diameter
- f) A breast tumor that is less than 5 mm away from the skin and chest wall, unless the tumor may be moved more than 5 mm with a bolus injection of saline or local anesthetic
- g) Patients unable to lie supine for up to 30 minutes

5 **Warnings**

Physicians are solely responsible for all clinical use of the Novilase Laser Therapy system and any results obtained with the device.

The Novilase Laser Therapy System should be used by physicians who have experience performing image-guided breast biopsy procedures or have experience with ultrasound image-guided techniques and have received appropriate training on the Novilase device.

When treating malignant breast tumors with the Novilase Laser Therapy system, MRI must be used to evaluate the absence of residual tumor at four to six weeks following the Novilase procedure. A breast biopsy may be performed to confirm malignancy if there is suspicion of a residual tumor. If necessary, retreatment may be performed with Novilase or with surgery.

6 **Precautions**

The Novilase Breast Laser Therapy Kit is intended for single use only. Do NOT re-sterilize or reuse. Possible risks associated with the reuse of a single-use device include but are not limited to mechanical failure, material degradation, and transmission of infectious agents.

- a) The Novilase laser and thermal probes are not compatible with magnetic resonance imaging for image guidance.
- b) Do not use open or damaged Novilase Breast ILT Kit pouches.
- c) Inspect devices upon removal of the protective packaging. Do not use a kinked fiber or damaged devices. Return damaged devices and packaging to the manufacturer or distributor.

- d) Use before the "Use By" date on the label.
- e) A spare treatment kit should be available in the procedure room before starting the Novilase procedure.
- f) Handle probes and the optical fiber as sharp objects.
- g) Local anesthetic may be used at the physician's discretion under the approved uses of this system. Use of more than 40 ml of local anesthetic and saline may result in suboptimal heating.
- h) Upon removal from the protective packaging, the laser fiber will tend to expand to full length. The user should be careful not to bend, kink, or damage the fiber. Be careful not touch the end faces of the fiber on any surface.
- i) Verification of focal destruction of the tumor (i.e., absence of residual tumor) may be obtained by MRI four weeks following the Novilase procedure. The subject's movement during a breast MRI exam should be limited. Ideally, the same MRI unit should be used for pre- and post-procedure imaging.
- j) The Novilase software requires a unique login and password to access and run. It runs on a MS Windows 10 or better operating system. No other software may be installed onto the system.

7 Residual Risks and undesirable side-effects

Predictable adverse events related to the Novilase device and procedure are similar to those experienced during a breast biopsy, including, but not limited to:

- a) Hemorrhage, hematoma, seroma, bruising, local swelling
- b) Infections
- c) Adverse reactions to anesthesia
- d) Pain, discomfort
- e) Minor skin scalding, blister
- f) Local redness of the skin
- g) Slight risk of seeding of tumor cells along the treatment channel.
- h) Microscopic tumor foci might remain after the procedure, and radiation and adjuvant therapy should be completed as prescribed
- i) Fat necrosis

Users and/or patients must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the country in which the user and/or patient is established.

8 Storage and handling conditions

Store the Novilase Breast ILT Kit in original packaging at room temperature. Following a procedure, dispose of probes, stylet, optical fiber, and all other sharp objects in approved sharp object disposal containers.

9 Novilase Procedure

9.1 Set up Saline Infusion Pump

- a) Refer to the Novilase LTS 2.1 User Manual and instructions for use for the respective saline pump.

9.2 Assemble Novilase Breast ILT Kit components before a patient enters the room

- a) Drape a sterile cover over Mayo stand or tray.
- b) Open a Novilase Laser Therapy Pack and remove the Novilase Breast ILT Kit and syringe. Open the sterile pouches and place contents onto a sterile tray. Remove probes from backer card. Place sterile scalpel on the tray.
- c) Enter Novilase Breast ILT Kit lot, serial, and control numbers from the kit label in Novilase system software to register the kit for use.
- d) Attach laser and thermal probes to their respective color-coded connector ports, using each port's keyhole as a guide.
- e) Remove protective cover from laser fiber connector. Connect laser fiber to the connector on the umbilical arm.
- f) Remove the stylet from the laser probe. Connect hemostasis valve to laser probe. Remove laser fiber from the protective tube. Insert fiber through the hemostasis valve until it protrudes past the tip of the cannula by 1 mm. Tighten the top valve to secure fiber in place.

- g) Press Test then Next.

System Test

Patient ID: MFG-005
Patient Name: Beatrice, Beatrice F

TREATMENT KIT

Lot Number: 0830
Control Number: 200004
Serial Number: 7F90B81F

DEVICE STATUS

Laser Probe Thermistor	OK
Thermal Probe Thermistor 1:	OK
Thermal Probe Thermistor 2:	OK
Thermal Probe Thermistor 3:	OK
Thermal Probe Thermistor 4:	OK
Thermal Probe Thermistor 5:	OK
Laser	OK

NEXT **CANCEL**

- h) If the test does not pass, follow the directions in the pop-up message.
- i) If the test passes, verify and acknowledge that all steps in the safety checklist pop-up shown below have been performed.

Verification and Acknowledgment

Verify and Acknowledge all check list item steps have been performed.
A new Novilase treatment kit is being used.
Disposable optical fiber has no kinks, tears, or bends tighter than a 6 inch diameter.
All disposable items are properly placed.
The umbilical has no kinks, tears, or bends tighter than a 6 inch diameter.
The infusion pump is set up and the saline tubing is acceptable.;
Bolus with sterile saline is finished.
All non-essential personnel have exited the treatment room.
All personnel in the room, including the patient, are wearing protective eyewear.
No error messages are displayed on the screen.
Press Ok to continue.

OK

- j) Remove the hemostasis valve by unscrewing the lower valve and removing fiber from the laser probe. Place this assembly in the sterile field for later use. Reinsert stylet.
- k) Fill syringe with sterile saline. Connect saline tube to hemostasis valve and the syringe. Place syringe onto the saline pump and set up saline pump per manufacturer's instructions. Bolus the saline line.
- l) Invite patient into room and position patient on table.

9.3 Use of Novilase disposables during the procedure

- a) Review pre-procedure imaging and locate the target tumor with ultrasound imaging. Any breast ultrasound with a 10 MHz or greater transducer may be used.

- b) Prep skin with Betadine or a suitable alternative. Inject local anesthesia.
- c) Capture pre-treatment images. Ensure that the tumor is 5 mm from the skin or chest wall or inject a bolus of local anesthetic to move the tumor to provide a 5 m separation.
- d) Nick the skin with a scalpel. Utilizing ultrasound image guidance, insert the laser probe with trocar stylet into the breast so that the tip is positioned in the center of the tumor or target ablation zone.
- e) Select the appropriate side of the probe guide based on the tumor's longest diameter and desired radius of the ablation zone.
 - i. Use the short side of probe guide for a tumor < 15 mm in longest diameter.
 - ii. Use the long side of the probe guide for a tumor 15 to 20 mm in longest diameter. The long side also may be chosen for smaller tumors to generate a larger ablative margin around the target tumor.

Snap the probe guide onto the laser probe flush with the shoulder on the probe.

- f) Nick the skin with a scalpel. Insert the thermal probe through the center of the probe guide until it is flush with the shoulder on the thermal probe, and parallel to the laser probe, as shown below.

- g) Remove the trocar stylet from the laser probe. Insert previously measured fiber through the cannula and tightly connect the hemostasis valve to the laser probe.
- h) Confirm the positioning of the probes with ultrasound imaging relative to the target tumor before initiating the procedure. Confirm that the laser probe is positioned in the target tumor or ablation zone center.
- i) Ensure everyone in the procedure room is wearing protective eyewear, including the patient. Ensure that a Laser in Use sign has been placed on door to ensure that no one enters procedure room during a procedure without protective eyewear.
- j) The Novilase procedure proceeds as described in the LTS 2.1 User Manual. The saline pump is operated according to the manufacturer's instructions for use.
- k) Select wattage to initiate ablation. Typically, start at 4 watts. In some cases, the entire power range of the Novilase system (i.e., up to 8 watts) may be required to reach target temperatures. Monitor the progress of the ablation by observing the temperature outputs

on the screen. Visualize the progress of the ablation on ultrasound. Target temperatures are achieved when the thermal probe sensors have reached 60°C or have exceeded 51°C for at least 3 minutes.

- l) Ensure target temperatures are achieved before pausing or ending the procedure and turning off the saline pump.
- m) Following the Novilase procedure, the probes and fiber are removed from the breast and discarded in a sharp's container.

10 Post-Novilase Procedure Instructions

These instructions are intended to serve as a guide. Physicians should use their judgment to recommend post-procedure instructions and follow-ups suited to their patients.


- a) Following the procedure, apply direct pressure to achieve hemostasis.
- b) Apply cold or ice packs to probe insertion site, wrap with a gauze bandage.
- c) Instruct patient to keep the cold pack on for 20 minutes and keep it off for 20 minutes. Repeat as necessary during the first day.
- d) Non-steroidal anti-inflammatory medications (NSAIDs) or acetaminophen may be recommended to control post-procedure discomfort or pain.
- e) Some bleeding at the site is to be expected. Instruct patients to call their physician if they still notice bleeding from the site after several hours have passed following the procedure or if they are in discomfort or notice excessive drainage from the site.
- f) Patients may resume activities of daily living shortly following the procedure but should avoid strenuous work or exercise for up to a week.


11 Expected Clinical Benefit

The Novilase Laser Therapy System is intended for use in ablating breast tumor tissue.

Compared to breast conserving surgery (BCS), the Novilase procedure provides an at least as good focal destruction of the breast tumor with minimal thermal damage to structures outside the targeted zone and a comparable safety profile but does not result in breast deformation (cosmetic outcome) and is performed percutaneously under local anesthesia within 30 minutes allowing for a fast recovery. The focal destruction of the tumor can be evaluated using MRI with a Negative Predictive Value not inferior to surgical margin pathology.

12 Symbols

Symbol	Meaning
	Manufacturer
	Authorized Representative in the European Community
	Authorized Representative in Switzerland
	Medical device per MDR
	Use-by date
	Serial number
	Lot number
	Control number, unique to each kit to prevent reuse
	Catalog number
	Consult instructions for use
	Do not reuse; single-use
	Do not re-sterilize
	Sterilization using irradiation
	Do not use if the package is damaged

Symbol	Meaning
	MRI unsafe. Steel probes are not compatible with MRI.
	Prescription only - device restricted to use by or on the order of a physician
	Single sterile barrier system with protective packaging inside
	Single sterile barrier system
	Laser in Use
	AC Input

Novilase® Brust-ILT-Kit

Modellnummer: 100010-01 und 100010-02

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	<p>Novian Health Inc. 222 W. Merchandise Mart Plaza Suite 570 Chicago, IL 60654 USA 312-266-7200</p>
---	--

E.C. Rep Ltd

5 Fitzwilliam Square East
Dublin 2, D02 R744, Ireland
Tel: +353 89 225 1951
Email: info@ecrep.ie

UK Responsible Person

Responsible Person Ltd.
Summit House, 4-5 Mitchell Street
Edinburgh, EH6 7BD Scotland
+44 131 618 1115

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